

OverCEM SA Overfibers Posts

Features Advantages

- Active Adhesion System*** Increases the adhesion with OverCEM SA
- Highly retentive surface** Increases the post/cement adhesion
- No silanization need** Time-saving and cost-effective
- Glass fibers** High flexure strength
- Balanced flexure E-modulus** To better stabilize the restoration
- Translucent** To accelerate the activation of the dual-cement
- Radiopaque** To be clearly distinguishable inside the canal
- Easy to remove**** Safely with absolute certainty and without tooth damages



OverCEM SA

Dual-Curing Self-Adhesive Resin Cement

- Self-Adhesive
- Innovative Formula
- Strong
- Dual-curing
- Auto-mixing
- T-Mixers tips
- Multi-purpose
- Luting cement




SIMPLE STRONG RELIABLE

www.overfibers.com



CODE	Description	OverCEM SA	T-Mixer Regular	T-Mixer Wide + Endo Tips
SKOCS	OverCEM SA Starter KIT	1 Translucent 9,5g syringe 1 Universal 9,5g syringe 1 Opaque 9,5g syringe	24 T-Mixer Regular	15 T-Mixer Wide + 15 Endo Tips
KOCTR	OverCEM SA Translucent KIT	3 Translucent 9,5g syringe	24 T-Mixer Regular	15 T-Mixer Wide + 15 Endo Tips
KOCUN	OverCEM SA Universal KIT	3 Universal 9,5g syringe	24 T-Mixer Regular	15 T-Mixer Wide + 15 Endo Tips
KOCOP	OverCEM SA Opaque KIT	3 Opaque 9,5g syringe	24 T-Mixer Regular	15 T-Mixer Wide + 15 Endo Tips
ROCTR	OverCEM SA Translucent	1 Translucent 9,5g syringe	8 T-Mixer Regular	5 T-Mixer Wide + 5 Endo Tips
ROCUN	OverCEM SA Universal	1 Universal 9,5g syringe	8 T-Mixer Regular	5 T-Mixer Wide + 5 Endo Tips
ROCOF	OverCEM SA Opaque	1 Opaque 9,5g syringe	8 T-Mixer Regular	5 T-Mixer Wide + 5 Endo Tips
TMREG	OverCEM SA T-Mixer Regular	25 T-Mixer Regular		
TMWE	OverCEM SA T-Mixer Wide + Endo Tip	15 T-Mixer Wide + 15 Endo Tips		

*Hi-Rem+ and Over Post+
**Hi-Rem+ and Hi-Rem Post



5 Posts for each size 0-1-2-3
1 Drill for each size 0-1-2-3

CODE	Description	OverCEM SA	T-Mixer Regular	T-Mixer Wide + 5 Endo Tips	Posts	Drills
SKOCH+	OverCEM SA - Hi-Rem+	1 Translucent 9,5g syringe	8 T-Mixer Regular	5 T-Mixer Wide + 5 Endo Tips	20 Posts	4 Drills
SKOCHP	OverCEM SA - Hi-Rem Prosthetic Post	1 Translucent 9,5g syringe	8 T-Mixer Regular	5 T-Mixer Wide + 5 Endo Tips	20 Posts	4 Drills
SKOCP+	OverCEM SA - Over Post+	1 Translucent 9,5g syringe	8 T-Mixer Regular	5 T-Mixer Wide + 5 Endo Tips	20 Posts	4 Drills
SKOCP+	OverCEM SA - Prosthetic Over Post	1 Translucent 9,5g syringe	8 T-Mixer Regular	5 T-Mixer Wide + 5 Endo Tips	20 Posts	4 Drills

Hi-Rem+

Hi-Rem Prosthetic Post

Over Post+

Prosthetic Over Post

Drill

CODE	Description	CODE	Description	CODE	Description	CODE	Description	Size Ø	CODE
RHP600	10 Hi-Rem+ 0	RPOP100	10 Hi-Rem Prosthetic Post 0	PP700	10 Over Post+ 0	POP100	10 Prosthetic Over Post 0	0.70 - 1.20	PD100
RHP601	10 Hi-Rem+ 1	RPOP101	10 Hi-Rem Prosthetic Post 1	PP701	10 Over Post+ 1	POP101	10 Prosthetic Over Post 1	0.80 - 1.40	PD101
RHP602	10 Hi-Rem+ 2	RPOP102	10 Hi-Rem Prosthetic Post 2	PP702	10 Over Post+ 2	POP102	10 Prosthetic Over Post 2	0.80 - 1.60	PD102
RHP603	10 Hi-Rem+ 3	RPOP103	10 Hi-Rem Prosthetic Post 3	PP703	10 Over Post+ 3	POP103	10 Prosthetic Over Post 3	0.90 - 1.80	PD103
						POP104	10 Prosthetic Over Post 4	1.10 - 2.10	PD104

OverCEM SA Starter KIT

OverCEM SA KIT

OverCEM SA Refill

T-Mixer Refill

OverCEM SA & Overfibers Posts

Overfibers Posts - Refill

Overfibers S.r.l.
Via E. Malatesta,7 - 40026 Imola (BO) - Italy
Phone +39 0542 52 153 - Fax +39 fax 0542 190 00 78
www.overfibers.com - info@overfibers.com



22 10 2020



OverCEM SA

Dual-Curing Self-Adhesive Resin Cement

Self-Adhesive Innovative Formula

Excellent bond strength stable over time

Strong

High Mechanical Properties

Dual-curing

Guarantees high polymerization efficiency


Auto-mixing

T-Mixer

-30% cement waste compared to conventional helix mixing tips

Multi-purpose Luting cement


Zirconia, Ceramics, Metals, Lithium disilicate, Posts, Composites



OverCEM SA Translucent OverCEM SA Universal OverCEM SA Opaque


Excellent bond strength to Zirconia and Lithium-disilicate ceramics

OverCEM SA reaches among most high bonding values of cements to Cubic-zirconia and Lithium-disilicate ceramics and optimal to Dentin.



Parisi C, Greischberger T, Baldissara P, Melilli D, Brusca G, Monaco C. Micro Shear Bond Strength of 10-MDP and 4-META Resin Cements. IADR General Session 2019 Vancouver, BC, Canada Abst. #1920

Outstanding Adhesion to Hi-Rem+ and Over Post+ fiber posts!



The cement/post interface showed high push-out shear forces: the adhesion between **Hi-Rem+** and **OverCEM SA** is increased thanks to the innovative **Active Adhesion System** contained in the post's resin matrix, and the extreme micro-roughness of its surface, greater than 5.0 microns (Ra).

OverCEM SA cement also showed a strong adhesion to the dentin, the micrograph shows the balanced adhesion between the two interfaces: cement/dentin and cement/post.

OverCEM SA

Dual-Curing Self-Adhesive Resin Cement

Features


Advantages

- T-Mixer the smallest on the market → Less cement waste in mixing tip
- Ideal consistency → More practical
- No etching - No priming - No bonding → Time and money saving
- Easy removal of excess → Faster use
- No post-operative sensitivity → No discomfort for the patient
- Strong and moisture-tolerant → Longer stability in time

High chromatic stability for outstanding aesthetics

OverCEM SA has the lowest absorption result

Before



After

OverCEM SA Universal Competitor 1 Competitor 2 Competitor 3




Even after 3 days at 37°C in concentrated coffee OverCEM SA color changing is negligible in comparison to other SA cements

Source Overfibers internal data

Indications for use

OverCEM SA Specific when you need

STRONG ADHESION



	Crowns	Bridges	Implant prosthodontics (Screwed metal-bonded)	Veneer crowns	Inlays/Onlays	Posts and dowels
Metal	●	●	●		●	●
Glass-ceramics	●	●	●	●	●	
Zirconia	●	●	●	●	●	●
Composites	●	●	●	●	●	●
Titanium Ceramic, Co-Cr Ceramic (implant prosthodontics)	●	●	●			